



ADJUSTABLE RESTING HAND SPLINTS

The resting hand splints are adjustable at the wrist to accommodate wrist contractures as well as hand/finger contractures. These splints are ideal for severe contractures.



**ADJUSTABLE GRIP HAND
ORTHOSIS SPLINT**



**ADJUSTABLE HAND THUMB
RESTING HAND SPLINT**

Comes with a below the wrist swivel that allows medial and lateral positions. The adjustable grip splint comes with 2 graduated finger rolls.



Brace moves left & right at the wrist.

see the catalog page 14 for more information...

Comes in 2 variations:

1. DEVIATION MODEL

corrects and controls radial or ulnar deviation. (Brace moves left and right at the wrist)

2. GONIOMETER MODEL

(Brace moves up and down at the wrist) can lock the wrist in any one of 7 possible positions. From 0 degrees to 180 degrees, the Goniometer Hand Thumb Resting Hand Splint can be graduated up to increase finger and wrist extension.

see the catalog page 15 for more information...

