

## Meet the Pneumatic Compression Device



Do you have patients suffering from Lymphedema or venous Insufficiency? Do they want a portable treatment that can be done at their own leisure?

This circulation pump promotes lymphatic and vascular flow by moving the fluid in a distal to proximal direction. The pump reduces swelling, promotes healing and improves range of motion in the affected limb(s).



# QUICK DELIVERY COVERED BY PATIENT'S INSURANCE

### **Gradient Sequential Compression Device**

HCPCS E0651

### **Device Features**

- Gradient Sequential Compression Device Most firm at the bottom (foot) and gradually reduces toward the top of the garment like a human body.
- Prevents blood and lymph fluid from reversing to leg or arm during inflation time.
- Gradient 4 Chamber Sequential Compression Therapy
- Pressure Adjustable from 20 mmHg ~ 100 mmHg
- Timer: 0 ~ 60 min & Continuous Mode

## **Gradient Sequential Compression Garments**

HCPCS E0667, E0668

### **Available Garments**

Full arm sleeve with shoulder and strap

Full arm sleeve w/o shoulder strap

Full leg sleeve (LG, XL, XXL)

Waist sleeve (LG, XL)

Medicare covers Pneumatic compression devices/Lymphedema Pumps if the following criteria is met.

- Patient has Chronic venous insufficiency (CVI) WITH venous stasis ulcers AND patient was treated for 6 months by physician with compression bandages or garments.
- Patient has one of the following: 1. Patient had radical cancer surgery or radiation therapy that interrupted the normal lymphatic fluid drainage. 2. Patient has malignant tumor with obstruction of the lymphatic drainage. 3. Patient has lymphedema since childhood or adolescence.

  If the above criteria is not met, Medicare will deny the device as not medically necessary